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*** Introduction**

Feline leukemia virus (FeLV) is a typical C-type oncovirus. It is the main cause of death among cats. It mainly affects cats under 5 years old and overcrowded. It spreads through licking, biting and ingestion of excretions and secretions. A 27 kilodalton protein moiety, p27, is a structural component of the inner viral core and is the major FeLV group-specific antigen. It can be found in great abundance in the cytoplasm of infected leukocytes and platelets and in soluble form in the plasma and serum of viremic cats. Positive detection of protein p27, the positive result, implicates that at least infected blood cells have shed the FeLV p27 in the circulating blood fraction.

Feline Immunodeficiency Virus (FIV) develops a slow infection, characterized by mainly affecting lymphocytes and macrophages, particularly CD4-marked lymphocytes, resembling the human immunodeficiency virus. It usually appears after 5 years of age due to its prolonged incubation period. It can also be transmitted through close and prolonged contact and by transplacental infection and through colostrum. Once the individual is infected it remains as a carrier for life. FIV is isolated in saliva, blood, tears and cerebrospinal fluid. From the infected cats two thirds do not develop the disease cause by this virus (unlike the HIV), in the remaining third the virus modifies the immune response, causing the immunodeficiency that eventually leads to secondary infections that may cause animal death.

*** Principles**

SensPERT™ FeLV Ag / FIV Ab test kit is for the simultaneous detection of FeLV p27 antigens and anti-FIV antibodies in whole blood, serum or plasma. The test results can appear on Control (C) and Test (T) lines where the principles of immunochromatography are used.

*** Characteristics**

- 1) One-step rapid test of FeLV Ag / FIV Ab
- 2) Rapid test results between 5 ~ 10 minutes
- 3) Expensive equipment not required
- 4) Easy storage and maintenance
- 5) High-purity and high-quality materials of the test kit increase its sensitivity and specificity.

*** Materials (10 tests/kit)**

- 1) Test 10 units
- 2) Diluent (buffer) 3 ml x 1 unit
- 3) Disposable dropper for sample collection 10 units

*** Composition**

Specimen well (S: for dropping), Test line (T), and Control line (C) are marked on the device. Inside of it, the strip is composed of sample pad, conjugate pad, nitrocellulose membrane (test paper) and absorbent pad

*** Effect**

Simultaneous detection of FeLV p27 antigens and anti-FIV antibodies in blood, serum or plasma

*** Uses**

- 1) **Specimen**
Feline whole blood, serum or plasma
- 2) **Test procedure**
 - a) When specimen and test kits are stored at cold circumstances (2~8℃), put them at room temperature for 15~30 minutes before use.
 - b) Take out a device from a pouch and place it on a horizontal surface.
 - c) Using a dropper as a pipette, obtain the specimen and dispense 2 drops (20 μl) of the fluid into the FeLV Ag specimen well (S) and 1 drop (10μl) of the fluid into the FIV Ab specimen well (S) of appropriately labeled device.
 - d) When the specimen is completely absorbed into each the specimen well (S), drop about 2 drops (80μl) of buffer, respectively.
 - e) Read the test results between 5 ~ 10 minutes.


◆ Invalid result after 10 min

3) Interpretation of the results

A purple band should appear on the control line regardless of the test result. The presence of another band on the test line determines the result.

Control Line (C): The line should always appear. If this line does not appear, the test should be considered invalid. This might be because of impure buffer or the lack of specimen. It should be tested again with new buffer.

Test Line (T): The presence of FeLV Ag and/or FIV Ab determines the test line.

2) FeLV Positive / FIV Negative

3) FeLV Negative / FIV Positive

4) FeLV / FIV Positive

5) Retesting

a) Both Test and Control lines do not appear.

b) Test line only appears.

FeLV Positive (+):

A positive result signifies the presence of circulating FeLV p27 in the blood fraction (serum, plasma, or whole blood) tested. Most but not all cats that test positive by SensPERT™ are actively excreting FeLV. It is thus important that positive FeLV tests be repeated within 2 to 3 months in order to determine whether the viremia is transient or persistent.

FeLV Negative (-):

A negative result indicates that no detectable p27 is present, but it does not exclude the possibility of viral incubation and is not an indication of immunity to FeLV. FeLV-negative infected with FeLV are not viremic and thus, do not shed infectious FeLV into their environment. Most latent infections appear to dissipate within a relatively short period of time (less than a year), perhaps because the infected cells differentiate to extinction or are scavenged by immunological processes.

FIV Positive (+):

A positive result implies the cat has the antibody to FIV. The fact of positive result does not necessarily mean that it is in the infected status. If the cat had been vaccinated 4~6 weeks ago, the possibility of interference cannot be eliminated. If the cat is under 6 months old, it could be maternal antibodies.

FIV Negative (-):

If the cat has no antibody to FIV, it would be better to plan of vaccination. If the cat is suspected to be infected but has negative result, it need to be tested repeatedly 2 or 3 times more. The time of the test could have been the period of seroconversion.

*** Further examinations**

This test is for primary screening only. Consult veterinarians for further necessary examinations to obtain clinical test results.

*** Precautions**

- 1) Use for feline *in-vitro* diagnostic purposes only.
- 2) Use within 10 minutes after opening the pouch because the test kit is very sensitive to moisture and its effect may diminish.
- 3) Be careful of not touching the result window.
- 4) Every specimen should be used with different droppers.
- 5) For testing, the buffer included should be used.
- 6) Do not use specimen showing hemolysis or being contaminated with microbes, which may cause false positive or false negative result.
- 7) Deal with specimen carefully. They can deliver unknown virus or infectious bacteria.
- 8) Use disposable gloves when you suspect the infection caused by specimen. And wash your hands later.
- 9) Dispose solid wastes after sterilizing them at 121℃ for over 1 hour.
- 10) Do not use the kit when its pouch is torn, sealing is not good in shape or expiration date is passed.

*** Storage method**

The test kits stored at 2~30℃ can be used for 24 months after manufacturing. Do not keep them in a refrigerator. However, if they are stored under cold circumstances, keep them at room temperature for 15 ~ 30 minutes before use.

*** Exchange**

The test kits are manufactured under strict quality control system. Nevertheless, if they are deteriorated during delivery, ask our distributors for exchange.

*** Liability**

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from use of this product.

Date issued : Sep. 08, 2016
 A106-01E-02

VETALL
 Manufactured by
VetAll Laboratories
 #B-305, Samsong Techno Valley, 140, Tongil-ro,
 Deogyang-gu, Goyang-si, Gyeonggi-do, 10594, Korea
 Tel : 82-2-22193458, Fax : 82-2-22193457


Instruction for FeLV Ag / FIV Ab test

Sample fill line (black line) : 1 drop (10 μl)


1 Sampling & loading


2 Absorption

3 Buffer addition

4 Reading

5 Discard