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(Feline Parvovirus Test Kit)

## Principles

SensPERT Feline Parvovirus Test Kit is designed to detect the antigens of feline parvovirus in feline feces. Two monoclonal antibodies in the kit specifically bind with different epitopes of the antigens. After being absorbed into the cellulose pad, the antigens of feline parvovirus move and bind with gold-colloid complex of monoclonal anti-feline parvovirus of the conjugate pad, forming Ab-Ag complex. This complex, then, forms Ab-Ag-Ab direct sandwich binding with the antibody of another monoclonal anti-feline parvovirus in the nitrocellulose membrane. The test results can appear on Control and Test lines where the principles of immunochromatography are applied.

## Characteristics

- 1) One-step rapid test of feline parvovirus antigens
- 2) Rapid test results between 5 ~ 10 minutes
- 3) Expensive equipment not required
- 4) Easy storage and maintenance
- 5) High-purity and high-quality materials of the test kit increases its sensitivity and specificity.

## Materials (10 tests/kit)

- 1) Test ..... 10 units
- 2) Diluent (buffer) ..... 1 ml x 10 units
- 3) Disposable dropper ..... 10 units
- 4) Swabs for sample collection ..... 10 units

## Composition

Specimen well (S : for dropping), Test line (T), and Control line (C) are marked on the device. Inside of it, the strip is composed of sample pad, conjugate pad, nitrocellulose membrane (test paper) and absorbent pad

## Effect

Detection of feline parvovirus antigens from feline feces

## Uses

- 1) **Specimen**  
Feline feces
- 2) **Storage of specimen**
  - a) Store specimen at 2~8°C for use within 24 hours or at -25°C for later use.
  - b) Keep the specimen at 22~25°C before use.
- 3) **Test procedure**
  - a) Collect feline feces with a swab and put the swab into a sample tube containing 1 ml of buffer.  
**Do not use sticky or a lump of feces.**
  - b) Swirl the buffer with the swab.
  - c) Drop 4 drops (100 µl) into the specimen well (S).
  - d) Read the test results within 5~10 minutes.

◆ Invalid result after 10 minutes

## 4) Interpretation of the results

A purple band should appear on the control line regardless of the test result. The presence of another band on the test line determines the result.

Control Line (C): The line should always appear regardless of the presence of the antigens of feline parvovirus. If this line does not appear, the test should be considered invalid. This might be because of impure

or the lack of specimen. It should be tested again with another kit.

Test Line (T): The presence of antigens of feline parvovirus determines the test line.

- 1) Negative: Control line only appears.

- 2) Positive: Both Test and Control lines appear.

- 3) Retesting

- a) Both Test and Control lines do not appear.

- b) Test line only appears.

## 5) Further examinations

This test is for primary screening only. Consult veterinarians for further necessary examinations to obtain clinical test results.

## Precautions

- 1) Use for feline in-vitro diagnostic purposes only.
- 2) Use within 10 minutes after opening the pouch because the test kit is very sensitive to moisture and its effect may diminish.
- 3) Be careful of not touching the result window.
- 4) Every specimen should be used with different droppers.
- 5) For testing, the buffer included should be used.
- 6) Do not use specimen being contaminated with microbes, which may cause false positive or false negative result.
- 7) Deal with specimen carefully. They can deliver unknown virus or infectious bacteria.
- 8) Use disposable gloves when you suspect the infection cause by specimen. And wash your hands later.
- 9) Dispose solid wastes after sterilizing them at 121°C for over 1 hour.
- 10) Do not use the kit when its pouch is torn, sealing is not good in shape or expiration date is passed.

## Storage method

The test kits stored at 2~30°C can be used for 24 months after manufacturing. Do not keep them in a refrigerator. However, if they are stored under cold circumstances, keep them at room temperature for 15 ~ 30 minutes before use.

## Exchange

The test kits are manufactured under strict quality control system. Nevertheless, if they are deteriorated during delivery, ask our distributors for exchange.

## Liability

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from use of this product.

Date issued : Feb 16, 2017  
A110-01E-01

**VETALL** Manufactured by  
**VetAll Laboratories**  
#B-305 Samsung Techno Valley, 140 Tongil-ro, Deogyang-gu,  
Goyang-si, Gyeonggi-do, 10594, Korea  
TEL. +82-2-2219-3459 FAX. +82-2-2219-3457


# Instruction for FPV Ag test

The amount of sample : 100 µl

## 1 Sampling

1. Do not use sticky or a lump of feces.
2. Do not collect too much feces.  
Maximum quantity : 30 mg  
(Fig : Swab with 30 mg of feces)

Swirl the buffer with the swab


## 2 Use the Supernatant layer

Wait until feces particles sink to the bottom


Supernatant layer for use  
Feces particle layer

## 3 Loading

- 4 drops (100 µl) diluted sample.  
**\*\*Do not overload the sample!!**

## 4 Reading

00:05:00~00:10:00

Read the result between 5~10 min

## 5 Discard

00:10:00

: Discard the used device  
**Invalid result after 10min**